
### DISPOSABLE INJECTION SYSTEM

Model	Tip Angle (°)	Tip Diameter (mm)
LC16	45	1.6
LC24	45	2.4
Cart4SS	45	1.6
CartM	45	1.6

STERILE EO

Rx only

CE  
2797

PI33 Rev 7

**DEVICE DESCRIPTION**

Lenstec cartridges are sterile, single use, disposable accessories, used to fold the intraocular lens prior to implantation. The devices possess a "twin flap" design with a small tip to allow the surgeon to fold and insert the intraocular lens through a small incision.

**INTENDED USE**

The Lenstec family of cartridges is intended to aid qualified surgeons in the implantation of a foldable intraocular lenses into the eye following cataract extraction. These products are required to successfully facilitate implantation through a small incision without damaging the intraocular lens.

**IMPORTANT NOTICE**

It is highly recommended that the surgeon adhere to the recommendations, contraindications and warnings outlined in these instructions.

**RISKS**

The potential risks of using the Lenstec Disposable Injection Systems are as follows: Capsular damage, Endophthalmitis, Toxic anterior segment syndrome (TASS), Uveitis, Iritis, Anterior capsule fibrosis

**WARNINGS**

- Cartridge injection systems manufactured by Lenstec are intended to be used only with IOLs manufactured by Lenstec. A **Compatibility Guide** is provided at the end of this document.
- Surgical instructions for use of the cartridge injection system with Lenstec's IOLs are provided in the IOL instructions for use.
- Do not use the cartridge after the expiration date shown on the outside package label.
- Proper procedure is the responsibility of the individual surgeon.
- Medical facilities utilizing these devices, and their accessories (if any), must ensure proper disposal as medical waste.

**A, C, D**

The titanium injector and lens loader are reusable (following decontamination and sterilization). The injector and lens loader is supplied unsterile and must be cleaned and sterilized prior to initial use.

**Manual cleaning prior to and after use:**

1. Rinse injector or lens loader with deionized water and then clean the device and its crevices using a soft surgical/laboratory brush moistened with alcohol\*. Care should be taken to avoid damage to the extreme tip.
2. Inspect the instrument under magnification (10x) to check for gross contamination, wear or damage. Do not use/reuse if wear or damage are apparent. Re-clean if gross contamination is visible.

**Sterilization method(s) prior to and after use:**

- Place unwrapped instruments onto trays and steam sterilize using gravity displacement sterilizer for a minimum 135°C for an exposure time of 3 minutes and minimum dry time of 0-1 minute.  
OR
- Place wrapped instruments onto trays and steam sterilize using gravity displacement sterilizer for a minimum 132°C for an exposure time of 15 minutes and dry time of 15-30 minutes.  
OR
- Place wrapped instruments onto trays and steam sterilize using gravity displacement sterilizer for a minimum 121°C for an exposure time of 30 minutes and dry time of 15-30 minutes.

**B**

- Re-use of the disposable injector is strictly prohibited, as it raises serious safety and efficacy concerns.
- LENSTEC does not provide cleaning/sterilization instructions. An improperly cleaned and/or sterilized disposable injector can cause significant damage to a patient's vision and health, due in part to cross contamination induced infection.
- Once removed from its original packaging, the disposable injector can lose traceability. In the event the disposable injector is re-used, it is unlikely the user will know the correct expiry date or serial number.
- LENSTEC cannot guarantee stability or proper function of the disposable injector in the event that an injector is re-used.

**QUALITATIVE & QUANTITATIVE DATA FOR THE CARTRIDGE INJECTION SYSTEMS**

Cartridge Injection System: Polypropylene cartridge. The devices have been tested and proven safe in accordance with ISO 10993-5, ISO 10993-7, ISO 10993-10, ISO 10993-11.

**EXPIRATION DATE**

The expiration date on the lens package is the sterility expiration date.

**RETURNS POLICY**

Contact your Lenstec representative regarding the return goods policy.

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Bibliography

- Chang, D. F., Mamalis, N., "Guidelines for the Cleaning and Sterilization of Intracocular Surgical Instruments", Ophthalmic Instrument Cleaning & Sterilization Guidelines (OICS) Task Force

COMPATIBILITY GUIDE

MODEL	POWERS (D)	ACCESSORIES
 Softec HP1™	+10.0D to +30.0D (LC16; Cart45S)	
 Softec I™ Softec HD™ Softec HDY™	-5.0D to +26.0D (LC16; Cart45S)  -5.0D to +36.0D (LC24)	
 Softec HDO™	+5.0D to +26.0D (LC16; Cart45S)  +5.0D to +36.0D (LC24)	<b>I9011S Reusable Injectors</b> LC16 (1.6mm) LC24 (2.4mm) Cart45S (1.6mm) CartM (1.6mm)  <b>Disposable Injectors</b> LC16 (1.6mm) LC24 (2.4mm) Cart45SI (1.6mm) CartMI (1.6mm)
 SBL-3™ SBL-2™	+10.0D to +26.0D (LC16; Cart45S)  +10.0D to +36.0D (LC24)	
 Tetraflex HD™	+5.0D to +26.0D (LC16; Cart45S)  +5.0D to +36.0D (LC24)	